



# PowerSoil<sup>®</sup> DNA Isolation Kit

Catalog No.	Quantity
12888-50	50 Preps
12888-100	100 Preps

## *Instruction Manual*

New protocol instruction: *Shake Solution C4 to mix before using to ensure consistent results.*

*Inhibitor Removal Technology<sup>®</sup> (IRT) is a registered trademark of MO BIO Laboratories, Inc. and is covered by US patent protection as well as international patents pending.*



**Please recycle**

Version: 08202010

Technical Information: Toll free 1-800-606-6246, or 1-760-929-9911 Email: [technical@mobio.com](mailto:technical@mobio.com) Website: [www.mobio.com](http://www.mobio.com)



## Table of Contents

Introduction .....	3
Protocol Overview.....	3
Flow Chart.....	5
Equipment Required .....	6
Kit Contents & Storage .....	6
Precautions & Warnings .....	6
Protocols:	
Experienced User Protocol.....	7
Detailed Protocol (Describes what is happening at each step) .....	8
Vacuum Manifold Protocol .....	11
Hints & Troubleshooting Guide.....	13
Contact Information.....	15
Other Quality Products Available. ....	16

## Introduction

The PowerSoil<sup>®</sup> DNA Isolation Kit is comprised of a novel and proprietary method for isolating genomic DNA from environmental samples utilizing our patented Inhibitor Removal Technology<sup>®</sup> (IRT). The kit is intended for use with environmental samples containing a high humic acid content including difficult soil types such as compost, sediment, and manure. Other more common soil types have also been used successfully with this kit. The isolated DNA has a high level of purity allowing for more successful PCR amplification of organisms from the sample. PCR analysis has been performed to detect a variety of organisms including bacteria (e.g. *Bacillus subtilis*, *Bacillus anthracis*), fungi (e.g. yeasts, molds), algae and Actinomycetes (e.g. *Streptomyces*).

## Protocol Overview

The PowerSoil<sup>®</sup> DNA Isolation Kit distinguishes itself from MO BIO's UltraClean<sup>®</sup> Soil DNA Isolation Kit with a humic substance/brown color removal procedure. This procedure is effective at removing PCR inhibitors from even the most difficult soil types. Environmental samples are added to a bead beating tube for rapid and thorough homogenization. Cell lysis occurs by mechanical and chemical methods. Total genomic DNA is captured on a silica membrane in a spin column format. DNA is then washed and eluted from the membrane. DNA is then ready for PCR analysis and other downstream applications.

## Bead Beating Options

The PowerSoil<sup>®</sup> DNA Isolation Kit does not require homogenization using a high velocity bead beater. However, if the microorganism of interest requires stronger homogenization than provided by a vortex, or if using a bead beater is desired, the PowerSoil<sup>®</sup> DNA Isolation Kit may be used in conjunction with the PowerLyzer<sup>™</sup> 24 homogenizer. MO BIO now offers the PowerLyzer<sup>™</sup> PowerSoil<sup>®</sup> DNA Isolation Kit (cat# 12855-50) with a Bead Tube suitable for high powered bead beating of soil. For more information about these products, or for references using the PowerSoil<sup>®</sup> DNA Isolation Kit with a FastPrep<sup>®</sup> instrument, please contact Technical Service at 1-800-606-6246 or [technical@mobio.com](mailto:technical@mobio.com).

Additional information can be found at [www.mobio.com/blog](http://www.mobio.com/blog) in the following articles:  
<http://www.mobio.com/blog/2009/11/08/molecular-biology-of-soil-an-introduction/>  
<http://www.mobio.com/blog/2010/01/17/molecular-biology-of-soil-dna-isolation-part-i/>

## Optimized for complete homogenization of any sample



**PowerLyzer<sup>™</sup> 24**  
**Bench Top Bead-Based Homogenizer**  
**Catalog#13155**  
**([www.mobio.com/powerlyzer](http://www.mobio.com/powerlyzer))**



## PowerLyzer™ 24 Bench Top Bead-Based Homogenizer

The PowerLyzer™ 24 Bench Top Bead-Based Homogenizer is a bead beating instrument uniquely designed for the most efficient and complete lysis and homogenization of any biological sample. In as little as 30 seconds, the PowerLyzer™ 24 homogenizer is capable of processing up to 24 samples in 2 ml tubes. With true "hands-free" operation, the downtime associated with manipulating samples through multiple cycles is eliminated. Even the toughest and most difficult samples such as pine needles, seeds, spores, fungal mats, and clay soils are easily and effectively lysed. For more information and protocols, call technical service.


## High Throughput Options

MO BIO offers a vacuum based protocol for faster processing without centrifugation for the DNA binding and column washing steps for Spin Filters. The MO BIO PowerVac™ Manifold allows for processing of up to 20 spin filter preps at a time using the PowerVac™ Mini Spin Filter Adapters. For additional high throughput options MO BIO offers the PowerSoil®-htp 96 Well Soil DNA Isolation Kit for processing up to 2 x 96 samples using a centrifuge capable of spinning two 96 Well Blocks stacked (13 cm x 8 cm x 5.5 cm) at 2500 x *g*. For 96 well homogenization of soil, MO BIO offers the 96 Well Plate Shaker and Plate Adapter Set (MO BIO Catalog# 11996 & 11999, respectively.)

**This kit is for research purposes only. Not for diagnostic use.**

Other Related Products	Catalog No.	Quantity
PowerMax® Soil DNA Isolation Kit	12988-10	10 preps
PowerSoil®-htp 96 Well Soil DNA Isolation Kit	12955-4	4 x 96 preps
	12955-12	12 x 96 preps
Ceramic Bead Tubes, 1.4 mm	13113-50	50 tubes
Glass Bead Tubes, 0.5 mm	13116-50	50 tubes
Glass Bead Tubes, 0.1mm	13118-50	50 tubes
PowerVac™ Manifold	11991	1 manifold
PowerVac™ Mini System	11992	1 unit + 20 adapters
PowerVac™ Mini Spin Filter Adapters	11992-10	10 adapters
	11992-20	20 adapters

## PowerSoil<sup>®</sup> DNA Isolation Kit





### Equipment Required

Microcentrifuge (10,000 x g)  
 Pipettors (50 µl - 500 µl)  
 Vortex-Genie<sup>®</sup> 2 Vortex (MO BIO Catalog# 13111-V or 13111-V-220)  
 Vortex Adapter (MO BIO Catalog # 13000-V1)

### Reagents Required but not Included

100% ethanol (for the PowerVac<sup>™</sup> Manifold protocol only)

### Kit Contents

Component	Kit Catalog # 12888-50		Kit Catalog # 12888-100	
	Catalog #	Amount	Catalog #	Amount
PowerBead Tubes (contain 750 µl solution)	12888-50-PBT	50	12888-100-PBT	100
PowerSoil <sup>®</sup> Solution C1	12888-50-1	3.3 ml	12888-100-1	6.6 ml
PowerSoil <sup>®</sup> Solution C2	12888-50-2	14 ml	12888-100-2	28 ml
PowerSoil <sup>®</sup> Solution C3	12888-50-3	11 ml	12888-100-3	22 ml
PowerSoil <sup>®</sup> Solution C4	12888-50-4	72 ml	12888-100-4	144 ml
PowerSoil <sup>®</sup> Solution C5	12888-50-5	30 ml	12888-100-5	2 x 30 ml
PowerSoil <sup>®</sup> Solution C6	12888-50-6	6 ml	12888-100-6	12 ml
PowerSoil <sup>®</sup> Spin Filters (units in 2 ml tubes)	12888-50-SF	50	12888-100-SF	100
PowerSoil <sup>®</sup> 2 ml Collection Tubes	12888-50-T	200	12888-100-T	400

### Kit Storage

Kit reagents and components should be stored at room temperature (15-30°C).

### Precautions

Please wear gloves when using this product. Avoid all skin contact with kit reagents. In case of contact, wash thoroughly with water. Do not ingest. See Material Safety Data Sheets for emergency procedures in case of accidental ingestion or contact. All MSDS information is available upon request (760-929-9911) or at [www.mobio.com](http://www.mobio.com). Reagents labeled flammable should be kept away from open flames and sparks.

**WARNING:** Solution C5 contains ethanol. It is flammable. Do not use bleach to clean the inside of the PowerVac<sup>™</sup> Manifold or to rinse the PowerVac<sup>™</sup> Mini Spin Filter Adapters when attached to the manifold.

**IMPORTANT NOTE FOR USE:** Make sure the 2 ml PowerBead Tubes rotate freely in your centrifuge without rubbing. Shake to mix Solution C4 before use.



## Experienced User Protocol

Please wear gloves at all times

1. To the **PowerBead Tubes** provided, 0.25 grams of soil sample.
2. Gently vortex to mix.
3. **Check Solution C1.** If **Solution C1** is precipitated, heat solution to 60°C until dissolved before use.
4. Add 60 µl of **Solution C1** and invert several times or vortex briefly.
5. Secure **PowerBead Tubes** horizontally using the MO BIO Vortex Adapter tube holder for the vortex (MO BIO Catalog# 13000-V1) or secure tubes horizontally on a flat-bed vortex pad with tape. Vortex at maximum speed for 10 minutes.  
**Note:** If you are using the 24 place Vortex Adapter for more than 12 preps, increase the vortex time by 5-10 minutes.
6. Make sure the PowerBead Tubes rotate freely in your centrifuge without rubbing. Centrifuge tubes at 10,000 x g for 30 seconds at room temperature. **CAUTION:** Be sure not to exceed 10,000 x g or tubes may break.
7. Transfer the supernatant to a clean **2 ml Collection Tube** (provided).  
**Note:** Expect between 400 to 500 µl of supernatant. Supernatant may still contain some soil particles.
8. Add 250 µl of **Solution C2** and vortex for 5 seconds. Incubate at 4°C for 5 minutes.
9. Centrifuge the tubes at room temperature for 1 minute at 10,000 x g.
10. Avoiding the pellet, transfer up to, but no more than, 600 µl of supernatant to a clean **2 ml Collection Tube** (provided).
11. Add 200 µl of **Solution C3** and vortex briefly. Incubate at 4°C for 5 minutes.
12. Centrifuge the tubes at room temperature for 1 minute at 10,000 x g.
13. Avoiding the pellet, transfer up to, but no more than, 750 µl of supernatant into a clean **2 ml Collection Tube** (provided).
14. Shake to mix Solution C4 before use. Add 1200 µl of **Solution C4** to the supernatant and vortex for 5 seconds.
15. Load approximately 675 µl onto a **Spin Filter** and centrifuge at 10,000 x g for 1 minute at room temperature. Discard the flow through and add an additional 675 µl of supernatant to the **Spin Filter** and centrifuge at 10,000 x g for 1 minute at room temperature. Load the remaining supernatant onto the **Spin Filter** and centrifuge at 10,000 x g for 1 minute at room temperature.  
**Note:** A total of three loads for each sample processed are required.
16. Add 500 µl of **Solution C5** and centrifuge at room temperature for 30 seconds at 10,000 x g.
17. Discard the flow through.
18. Centrifuge again at room temperature for 1 minute at 10,000 x g.
19. Carefully place spin filter in a clean **2 ml Collection Tube** (provided). Avoid splashing any **Solution C5** onto the **Spin Filter**.
20. Add 100 µl of **Solution C6** to the center of the white filter membrane. Alternatively, sterile DNA-Free PCR Grade Water may be used for elution from the silica Spin Filter membrane at this step (MO BIO Catalog# 17000-10).
21. Centrifuge at room temperature for 30 seconds at 10,000 x g.
22. Discard the **Spin Filter**. The DNA in the tube is now ready for any downstream application. No further steps are required.

We recommend storing DNA frozen (-20° to -80°C). **Solution C6** contains no EDTA. To concentrate the DNA see the Hints & Troubleshooting Guide.

Thank you for choosing the PowerSoil® DNA Isolation Kit.

## Detailed Protocol (Describes what is happening at each step)

Please wear gloves at all times

1. To the **PowerBead Tubes** provided, add 0.25 grams of soil sample.

*What's happening: After your sample has been loaded into the PowerBead Tube, the next step is a homogenization and lysis procedure. The PowerBead Tube contains a buffer that will (a) help disperse the soil particles, (b) begin to dissolve humic acids and (c) protect nucleic acids from degradation.*

2. Gently vortex to mix.

*What's happening: Gentle vortexing mixes the components in the PowerBead Tube and begins to disperse the sample in the PowerBead Solution.*

3. **Check Solution C1.** If **Solution C1** is precipitated, heat solution to 60°C until the precipitate has dissolved before use.

*What's happening: Solution C1 contains SDS and other disruption agents required for complete cell lysis. In addition to aiding in cell lysis, SDS is an anionic detergent that breaks down fatty acids and lipids associated with the cell membrane of several organisms. If it gets cold, it will form a white precipitate in the bottle. Heating to 60°C will dissolve the SDS and will not harm the SDS or the other disruption agents. Solution C1 can be used while it is still warm.*

4. Add 60 µl of **Solution C1** and invert several times or vortex briefly.

5. Secure **PowerBead Tubes horizontally** using the MO BIO Vortex Adapter tube holder for the vortex (MO BIO Catalog# 13000-V1) or secure tubes horizontally on a flat-bed vortex pad with tape. Vortex at maximum speed for 10 minutes. **Note:** If you are using the 24 place Vortex Adapter for more than 12 preps, increase the vortex time by 5-10 minutes.

**Note:** *The vortexing step is critical for complete homogenization and cell lysis. Cells are lysed by a combination of chemical agents from steps 1-4 and mechanical shaking introduced at this step. By randomly shaking the beads in the presence of disruption agents, collision of the beads with microbial cells will cause the cells to break open.*

*What's happening: The MO BIO Vortex Adapter is designed to be a simple platform to facilitate keeping the tubes tightly attached to the vortex. It should be noted that although you can attach tubes with tape, often the tape becomes loose and not all tubes will shake evenly or efficiently. This may lead to inconsistent results or lower yields. Therefore, the use of the MO BIO Vortex Adapter is a highly recommended and cost effective way to obtain maximum DNA yields.*

6. Make sure the **PowerBead Tubes** rotate freely in your centrifuge without rubbing. Centrifuge tubes at 10,000 x g for 30 seconds at room temperature. **CAUTION:** Be sure not to exceed 10,000 x g or tubes may break.
7. Transfer the supernatant to a clean **2 ml Collection Tube** (provided).

**Note:** *Expect between 400 to 500 µl of supernatant at this step. The exact recovered volume depends on the absorbency of your starting material and is not critical for the procedure to be effective. The supernatant may be dark in appearance and still contain some soil particles. The presence of carry over soil or a dark color in the mixture is expected in many soil types at this step. Subsequent steps in the protocol will remove both carry over soil and coloration of the mixture.*





8. Add 250  $\mu$ l of **Solution C2** and vortex for 5 seconds. Incubate at 4°C for 5 minutes.

*What's happening: Solution C2 is patented Inhibitor Removal Technology<sup>®</sup> (IRT). It contains a reagent to precipitate non-DNA organic and inorganic material including humic substances, cell debris, and proteins. It is important to remove contaminating organic and inorganic matter that may reduce DNA purity and inhibit downstream DNA applications.*

9. Centrifuge the tubes at room temperature for 1 minute at 10,000 x g.

10. Avoiding the pellet, transfer up to 600  $\mu$ l of supernatant to a clean **2 ml Collection Tube** (provided).

*What's happening: The pellet at this point contains non-DNA organic and inorganic material including humic acid, cell debris, and proteins. For the best DNA yields, and quality, avoid transferring any of the pellet.*

11. Add 200  $\mu$ l of **Solution C3** and vortex briefly. Incubate at 4°C for 5 minutes.

*What's happening: Solution C3 is patented Inhibitor Removal Technology<sup>®</sup> (IRT) and is a second reagent to precipitate additional non-DNA organic and inorganic material including humic acid, cell debris, and proteins. It is important to remove contaminating organic and inorganic matter that may reduce DNA purity and inhibit downstream DNA applications.*

12. Centrifuge the tubes at room temperature for 1 minute at 10,000 x g.

13. Transfer up to 750  $\mu$ l of supernatant to a clean **2 ml Collection Tube** (provided).

*What's happening: The pellet at this point contains additional non-DNA organic and inorganic material including humic acid, cell debris, and proteins. For the best DNA yields, and quality, avoid transferring any of the pellet.*

14. Shake to mix Solution C4 before use. Add 1.2 ml of **Solution C4** to the supernatant (be careful solution doesn't exceed rim of tube) and vortex for 5 seconds.

*What's happening: Solution C4 is a high concentration salt solution. Since DNA binds tightly to silica at high salt concentrations, this will adjust the DNA solution salt concentrations to allow binding of DNA, but not non-DNA organic and inorganic material that may still be present at low levels, to the Spin Filters.*

15. Load approximately 675  $\mu$ l onto a **Spin Filter** and centrifuge at 10,000 x g for 1 minute at room temperature. Discard the flow through and add an additional 675  $\mu$ l of supernatant to the **Spin Filter** and centrifuge at 10,000 x g for 1 minute at room temperature. Load the remaining supernatant onto the **Spin Filter** and centrifuge at 10,000 x g for 1 minute at room temperature.

**Note:** A total of three loads for each sample processed are required.

*What's happening: DNA is selectively bound to the silica membrane in the Spin Filter device in the high salt solution. Contaminants pass through the filter membrane, leaving only DNA bound to the membrane.*

16. Add 500  $\mu$ l of **Solution C5** and centrifuge at room temperature for 30 seconds at 10,000 x g.

*What's happening: Solution C5 is an ethanol based wash solution used to further clean the DNA that is bound to the silica filter membrane in the Spin Filter. This wash solution removes residual salt, humic acid, and other contaminants while allowing the DNA to stay bound to the silica membrane.*



17. Discard the flow through from the **2 ml Collection Tube**.

*What's happening: This flow through fraction is just non-DNA organic and inorganic waste removed from the silica Spin Filter membrane by the ethanol wash solution.*

18. Centrifuge at room temperature for 1 minute at 10,000 x g.

*What's happening: This second spin removes residual Solution C5 (ethanol wash solution). It is critical to remove all traces of wash solution because the ethanol in Solution C5 can interfere with many downstream DNA applications such as PCR, restriction digests, and gel electrophoresis.*

19. Carefully place Spin Filter in a clean **2 ml Collection Tube** (provided). Avoid splashing any **Solution C5** onto the **Spin Filter**.

**Note:** *It is important to avoid any traces of the ethanol based wash solution.*

20. Add 100  $\mu$ l of **Solution C6** to the center of the white filter membrane.

**Note:** *Placing the Solution C6 (sterile elution buffer) in the center of the small white membrane will make sure the entire membrane is wetted. This will result in a more efficient and complete release of the DNA from the silica Spin Filter membrane. As Solution C6 (elution buffer) passes through the silica membrane, DNA that was bound in the presence of high salt is selectively released by Solution C6 (10 mM Tris) which lacks salt.*

*Alternatively, sterile DNA-Free PCR Grade Water may be used for DNA elution from the silica Spin Filter membrane at this step (MO BIO Catalog# 17000-10). Solution C6 contains no EDTA. If DNA degradation is a concern, Sterile TE may also be used instead of Solution C6 for elution of DNA from the Spin Filter.*

21. Centrifuge at room temperature for 30 seconds at 10,000 x g.

22. Discard the **Spin Filter**. The DNA in the tube is now ready for any downstream application. No further steps are required.

We recommend storing DNA frozen (-20° to -80°C). **Solution C6** does not contain any EDTA. To concentrate DNA see the Hints & Troubleshooting Guide.

**Thank you for choosing the PowerSoil<sup>®</sup> DNA Isolation Kit.**



## Vacuum Protocol using the PowerVac™ Manifold

Please wear gloves at all times

For each sample lysate, use one Spin Filter column. Keep the Spin Filter in the attached 2 ml Collection Tube and continue using the Collection Tube as a Spin Filter holder until needed for the Vacuum Manifold Protocol. Label each Collection Tube top and Spin Filter column to maintain sample identity. If the Spin Filter becomes clogged during the vacuum procedure, you can switch to the procedure for purification of the DNA by centrifugation.

You will need to provide 100% ethanol for step 4 of this protocol

1. For each prep, attach one aluminum **PowerVac™ Mini Spin Filter Adapter** (MO BIO Catalog# 11992-10 or 11992-20) into the Luer-Lok® fitting of one port in the manifold. Gently press a Spin Filter column into the PowerVac™ Mini Spin Filter Adapter until snugly in place. Ensure that all unused ports of the vacuum manifold are closed.  
**Note:** Aluminum PowerVac™ Mini Spin Filter Adapters are reusable.
2. Transfer 650 µl of prepared sample lysate (from step 14) to the **Spin Filter column**.
3. Turn on the vacuum source and open the stopcock of the port. Hold the tube in place when opening the stopcock to keep the spin filter steady. Allow the lysate to pass through the **Spin Filter column**. After the lysate has passed through the column completely, load again with the next 650 µl of lysate. Continue until all of the lysate has been loaded onto the **Spin Filter column**. Close the one-way Luer-Lok® stopcock of that port.  
**Note:** If Spin Filter Columns are filtering slowly, close the ports to samples that have completed filtering to increase the pressure to the other columns.
4. Load 800 µl of 100% ethanol into the Spin Filter so that it completely fills the column. Open the stopcock while holding the column steady. Allow the ethanol to pass through the column completely. Close the stopcock.
5. Add 500 µl of **Solution C5** to each Spin Filter. Open the Luer-Lok® stopcock and apply a vacuum until **Solution C5** has passed through the Spin Filter completely. Continue to pull a vacuum for another minute to dry the membrane. Close each port.
6. Turn off the vacuum source and open an unused port to vent the manifold. If all 20 ports are in use, break the vacuum at the source. Make certain that all vacuum pressure is released before performing the next step. It is important to turn off the vacuum at the source to prevent backflow into the columns.
7. Remove the **Spin Filter column** and place in the original labeled **2 ml Collection Tube**. Place into the centrifuge and spin at 13,000 × *g* for 1 minute to completely dry the membrane.
8. Transfer the **Spin Filter column** to a new **2 ml Collection Tube** and add 100 µl of **Solution C6** to the center of the white filter membrane. Alternatively, sterile DNA-Free PCR Grade Water may be used for elution from the silica **Spin Filter** membrane at this step (MO BIO Catalog # 17000-10).
9. Centrifuge at room temperature for 30 seconds at 10,000 × *g*.



10. Discard the **Spin Filter column**. The DNA in the tube is now ready for any downstream application. No further steps are required.

We recommend storing DNA frozen (-20° to -80°C). **Solution C6** contains no EDTA. To concentrate the DNA see the Hints & Troubleshooting Guide.

**Thank you for choosing the PowerSoil<sup>®</sup> DNA Isolation Kit.**



## Hints & Troubleshooting Guide

### **Amount of Soil to Process**

This kit is designed to process 0.25 grams of soil. For inquiries regarding the use of larger sample amounts, please contact technical support for suggestions. For wet soils, see information under "Wet Soil Sample" below.

### **Wet Soil Sample**

If soil sample is high in water content, remove contents from PowerBead Tube (beads and solution) and transfer into another sterile microcentrifuge tube (not provided). Add soil sample to PowerBead Tube and centrifuge at room temperature for 30 seconds at 10,000 x g. Remove as much liquid as possible with a pipet tip. Add beads and bead solution back to PowerBead Tube and follow protocol starting at step 2.

### **If DNA Does Not Amplify**

- Make sure to check DNA yields by gel electrophoresis or spectrophotometer reading. An excess amount of DNA will inhibit a PCR reaction.
- Diluting the template DNA should not be necessary with DNA isolated with the PowerSoil<sup>®</sup> DNA Isolation Kit; however, it should still be attempted.
- If DNA will still not amplify after trying the steps above, then PCR optimization (changing reaction conditions and primer choice) may be needed.

### **Eluted DNA Sample Is Brown**

We have not observed any coloration in DNA isolated using the PowerSoil<sup>®</sup> DNA Isolation Kit. If you observe coloration in your samples, please contact technical support for suggestions.

### **Alternative Lysis Methods**

- After adding Solution C1, vortex 3-4 seconds, then heat to 70°C for 5 minutes. Vortex 3-4 seconds. Heat another 5 minutes. Vortex 3-4 seconds. This alternative procedure will reduce shearing but may also reduce yield.
- If cells are difficult to lyse, a 10 minute incubation at 70°C, after adding Solution C1, can be performed. Follow by continuing with protocol step 5.

### **Concentrating the DNA**

The final volume of eluted DNA will be 100 µl. The DNA may be concentrated by adding 4 µl of 5 M NaCl and inverting 3-5 times to mix. Next, add 200 µl of 100% cold ethanol and invert 3-5 times to mix. Centrifuge at 10,000 x g for 5 minutes at room temperature. Decant all liquid. Remove residual ethanol in a speed vac, dessicator, or air dry. Resuspend precipitated DNA in sterile water or sterile 10 mM Tris.

### **DNA Floats Out of Well When Loaded on a Gel**

This usually occurs because residual Solution C5 remains in the final sample. Prevent this by being careful in step 19 not to transfer liquid onto the bottom of the spin filter basket. Ethanol precipitation (described in "Concentrating the DNA") is the best way to remove residual Solution C5.

### **Storing DNA**

DNA is eluted in Solution C6 (10 mM Tris) and must be stored at -20° to -80°C to prevent degradation. DNA can be eluted in TE without loss, but the EDTA may inhibit downstream reactions such as PCR and automated sequencing. DNA may also be eluted with sterile DNA-Free PCR Grade Water (MO BIO Catalog# 17000-10).



## Hints & Troubleshooting Guide cont.

### ***Cleaning of the PowerVac™ Mini Spin Filter Adapters***

It is recommended to rinse the PowerVac™ Mini Spin Filter Adapters promptly after use to avoid salt build up. To clean the PowerVac™ Mini Spin Filter Adapters, rinse each adapter with DI water followed by 70% ethanol and flush into the manifold base. Alternatively, remove the adapters and wash in laboratory detergent and DI water. PowerVac™ Mini Spin Filter Adapters may be autoclaved.

**Do not use bleach to clean the PowerVac™ Mini Spin Filter Adapters while attached to the PowerVac™ Manifold. Bleach should never be mixed with solutions containing guanidine and should not be used to clean the PowerVac™ Manifold. For more information on cleaning the PowerVac™ Manifold, please refer to the PowerVac™ Manifold manual.**



## Contact Information

### Technical Support:

Phone MO BIO Laboratories, Inc. Toll Free 800-606-6246, or 760-929-9911

Email: [technical@mobio.com](mailto:technical@mobio.com)

Fax: 760-929-0109

Mail: MO BIO Laboratories, Inc, 2746 Loker Ave West, Carlsbad, CA 92010

### Ordering Information:

Direct: Phone MO BIO Laboratories, Inc. Toll Free 800-606-6246, or 760-929-9911

Email: [orders@mobio.com](mailto:orders@mobio.com)

Fax: 760-929-0109

Mail: MO BIO Laboratories, Inc, 2746 Loker Ave West, Carlsbad, CA 92010

For the distributor nearest you, visit our website at [www.mobio.com/distributors](http://www.mobio.com/distributors)



## Other Quality Products Available from MO BIO Laboratories, Inc.

For more product and detailed information go to [www.mobio.com/catalog-request](http://www.mobio.com/catalog-request) to request a catalog.

DNA Purification and Gel Extraction	Catalog No.	Quantity
PowerClean® DNA Clean-Up Kit	12877-50	50 preps
UltraClean® 15 DNA Purification Kit	12100-300	300 preps
UltraClean® PCR Clean-Up Kit	12500-50 12500-100 12500-250	50 preps 100 preps 250 preps
UltraClean®-htp 96 Well PCR Clean-Up Kit	12596-4 12596-12	4 x 96 preps 12 x 96 preps
UltraClean® GelSpin® DNA Extraction Kit	12400-50 12400-100 12400-250	50 preps 100 preps 250 preps
Plasmid DNA Isolation	Catalog No.	Quantity
UltraClean® 6 Minute Mini Plasmid Prep Kit	12300-50 12300-100 12300-250	50 preps 100 preps 250 preps
UltraClean® Standard Mini Plasmid Prep Kit	12301-50 12301-100 12301-250	50 preps 100 preps 250 preps
UltraClean®-htp 96 Well Plasmid Prep Kit	12396-4 12396-12	4 x 96 preps 12 x 96 preps
UltraClean® Midi Plasmid Prep Kit	12700-20 12700-50	20 preps 50 preps
UltraClean® Maxi Plasmid Prep Kit	12600-10 12600-20	10 preps 20 preps
UltraClean® Endotoxin-Free Mini Plasmid Prep Kit	12311-100 12311-250	100 preps 250 preps
UltraClean® Endotoxin-Free Midi Plasmid Prep Kit	12711-10	10 preps
UltraClean® Endotoxin-Free Maxi Plasmid Prep Kit	12611-10	10 preps
UltraClean® Endotoxin Removal Kit	12615	1 kit
UltraClean® Endotoxin-Free Ethanol Precipitation Kit	12616	1 kit
UltraClean® Endotoxin Removal Reagent	12625-25	25 ml
Endotoxin-Free Sodium Chloride	12626-15	15 ml
Endotoxin-Free Centrifuge Tubes	12617-100 12618-50 12619-25	100 each/2 ml tubes 50 each/15 ml tubes 25 each/50 ml tubes
RNA Isolation	Catalog No.	Quantity
PowerLyzer™ UltraClean® Tissue & Cells RNA Isolation Kit	15055-50	50 preps
PowerLyzer™ UltraClean® Plant RNA Isolation Kit	13355-50	50 preps
PowerBiofilm™ RNA Isolation Kit	25000-50	50 preps
LifeGuard™ Soil Stabilization Solution	12868-10 12868-100 12868-1000 12868-7500	10 ml 100 ml 1 L 7.5 L
On-Spin Column DNase I Kit (RNase-Free)	15100-50	50 preps
BiOstic® Stabilized Blood RNA Isolation Kit	12231-20 12231-50 12231-100	20 preps 50 preps 100 preps
BiOstic® Blood Total RNA Isolation Kit	12230-20 12230-50	20 preps 50 preps

RNA Isolation ... Continued	Catalog No.	Quantity
RNA PowerSoil® DNA Elution Accessory Kit	12867-25	25 preps
RNA PowerSoil® Total RNA Isolation Kit	12866-25	25 preps
UltraClean® Microbial RNA Isolation Kit	15800-50 15800-250	50 preps 250 preps
UltraClean® Tissue & Cells RNA Isolation Kit	15000-50 15000-250	50 preps 250 preps
UltraClean® Plant RNA Isolation Kit	13300-20 13300-50	20 preps 50 preps
Genomic DNA Isolation	Catalog No.	Quantity
PowerLyzer™ PowerSoil® DNA Isolation Kit	12855-50	50 preps
PowerLyzer™ UltraClean® Microbial DNA Isolation Kit	12255-50	50 preps
PowerBiofilm™ DNA Isolation Kit	24000-50	50 preps
PowerFood™ Microbial DNA Isolation Kit	21000-50 21000-100	50 preps 100 preps
BiOstic® Bacteremia DNA Isolation Kit	12240-50	50 preps
BiOstic® FFPE Tissue DNA Isolation Kit	12250-50	50 preps
BiOstic® Paraffin Removal Reagent	12251-50	2 x 25 ml
PowerMax® Soil DNA Isolation Kit	12988-10	10 preps
PowerSoil® DNA Isolation Kit	12888-50 12888-100	50 preps 100 preps
PowerSoil®-htp 96 Well Soil DNA Isolation Kit	12955-4 12955-12	4 x 96 preps 12 x 96 preps
UltraClean® Soil DNA Isolation Kit	12800-50 12800-100	50 preps 100 preps
UltraClean®-htp 96 Well Soil DNA Isolation Kit	12896-4 12896-12	4 x 96 preps 12 x 96 preps
UltraClean® Mega Soil DNA Isolation Kit	12900-10	10 preps
PowerClean® DNA Clean-Up Kit	12877-50	50 preps
UltraClean® Fecal DNA Isolation Kit	12811-50 12811-100	50 preps 100 preps
PowerMicrobial® Midi DNA Isolation Kit	12225-25	25 preps
PowerMicrobial® Maxi DNA Isolation Kit	12226-25	25 preps
UltraClean® Microbial DNA Isolation Kit	12224-50 12224-250	50 preps 250 preps
UltraClean®-htp 96 Well Microbial DNA Isolation Kit	10196-4 10196-12	4 x 96 preps 12 x 96 preps
PowerPlant® DNA Isolation Kit	13200-50 13200-100	50 preps 100 preps
UltraClean® Plant DNA Isolation Kit	13000-50 13000-250	50 preps 250 preps





## Other Quality Products Available from MO BIO Laboratories, Inc.

For more product and detailed information go to [www.mobio.com/catalog-request](http://www.mobio.com/catalog-request) to request a catalog.

Genomic DNA Isolation ...Continued	Catalog No.	Quantity	Other Reagents and Lab Accessories	Catalog No.	Quantity
UltraClean®-htp 96 Well Plant DNA Isolation Kit	13096-4 13096-12	4 x 96 preps 12 x 96 preps	LB Agar (Lennox) Powder Growth Media, pH 7	12109-05 12109-1 12109-5	500 g 1 kg 5 kg
UltraClean® Tissue & Cells DNA Isolation Kit	12334-50 12334-250	50 preps 250 preps	Soybean-Casein Digest Medium (TSB), USP	12114-05 12114-1 12114-5	500 g 1 kg 5 kg
UltraClean®-htp 96 Well Tissue DNA Isolation Kit	12996-4 12996-12	4 x 96 preps 12 x 96 preps	Soybean-Casein Digest Agar Medium (TSA), USP	12115-05 12115-1 12115-5	500 g 1 kg 5 kg
UltraClean® Blood DNA Isolation Kit (Non-Spin)	12000-100	100 preps	Yeast Extract	12110-05 12110-1 12110-5	500 g 1 kg 5 kg
UltraClean® Blood DNA Isolation Kit (Processes 1,000 ml of Blood)	12000-1000	1 kit	Tryptone	12111-05 12111-1 12111-5	500 g 1 kg 5 kg
UltraClean® Blood DNA Isolation Kit Plus RNase (Processes 1,000 ml of Blood)	12002-1000	1 kit	Agar, Bacteriological Grade	12112-05 12112-1 12112-5	500 g 1 kg 5 kg
UltraClean® BloodSpin® DNA Isolation Kit	12200-50 12200-250	50 preps 250 preps	20 bp DNA Ladder	17020-40	40 µg
UltraClean®-htp 96 Well BloodSpin® DNA Isolation Kit	12296-4 12296-12	4 x 96 preps 12 x 96 preps	100 bp DNA Ladder	17100-40	40 µg
UltraClean® Forensic DNA Isolation Kit	14000-10 14000-20	10 isolations 20 isolations	1 kb DNA Ladder	17200-100	100 µg
PowerWater® DNA Isolation Kit	14900-50-NF 14900-50-22 14900-50-45  14900-100-NF 14900-100-22 14900-100-45	50 preps (No filters) (0.22 µm) 100 preps (0.45 µm)  100 preps (No filters) (0.22 µm) 100 preps (0.45 µm)	UltraClean® Agarose, Molecular Biology Grade	15003-50 15003-100 15003-500 15003-1000	50 g 100 g 500 g 1 kg
RapidWater™ DNA Isolation Kit	14810-50-NF 14810-50-22 14810-50-45  14810-100-NF 14810-100-22 14810-100-45	50 preps (No filters) (0.22 µm) 100 preps (0.45 µm)  100 preps (No filters) (0.22 µm) 100 preps (0.45 µm)	UltraClean® MS-8 Agarose	15515-50 15515-100 15515-500	50 g 100 g 500 g
UltraClean® Water DNA Isolation Kit (0.45µm filters)	14800-10 14800-25	10 preps 25 preps	UltraClean® Forensic Agarose	15505-50 15505-100 15505-500	50 g 100 g 500 g
UltraClean® Water DNA Isolation Kit (0.22 µm filters)	14880-10 14880-25	10 preps 25 preps	UltraClean® Low Melt Agarose	15005-50 15005-100 15005-500	50 g 100 g 500 g
UltraClean® Water DNA Isolation Kit (No filters)	14800-10-NF 14800-25-NF	10 preps 25 preps	UltraClean® Low Melt Sieve Agarose	15004-50 15004-100 15004-500	50 g 100 g 500 g
<b>Microbiological Culture Media</b>	<b>Catalog No.</b>	<b>Quantity</b>	Ethidium Bromide Solution	15006-1 15006-10	1 ml 10 ml
TB DRY® Powder Growth Media	12105-05 12105-1 12105-5	500 g 1 kg 5 kg	Ethidium Bromide Destaining Tea Bags	15007-25	25 bags
LB Broth Powder Growth Media, pH 7	12106-05 12106-1 12106-5	500 g 1 kg 5 kg	Bromophenol Blue Gel Loading Buffer	15008-1 15008-5	1 ml 5 x 1 ml
LB Agar Powder Growth Media, pH 7	12107-05 12107-1 12107-5	500 g 1 kg 5 kg	Bromophenol Blue/Xylene Cyanol Gel Loading Buffer	15009-1 15009-5	1 ml 5 x 1 ml
LB Broth (Lennox) Powder Growth Media, pH 7	12108-05 12108-1 12108-5	500 g 1 kg 5 kg	TAE Buffer, 50X (Tris-acetate-EDTA)	15001-100 15001-500 15001-1000	100 ml 500 ml 1 liter



## Other Quality Products Available from MO BIO Laboratories, Inc.

For more product and detailed information go to [www.mobio.com/catalog-request](http://www.mobio.com/catalog-request) to request a catalog.

Other Reagents and Lab Accessories... Continued	Catalog No.	Quantity	Instrumentation and Accessories	Catalog No.	Quantity
TBE Buffer, 10X (Tris-borate-EDTA)	15002-100 15002-500 15002-1000	100 ml 500 ml 1 liter	PowerLyzer™ 24 Bench Top Bead-Based Homogenizer (110/220V)	13155	1 unit
RNase-Free Gloves	1555-XS 1555-S 1555-M 1555-L	bag of 100 bag of 100 bag of 100 bag of 100	PowerLyzer™ Tube Holder	13156	1 unit
UltraClean® Lab Cleaner	12095-250  12095-500  12095-1000	250 ml squeeze bottle 500 ml spray bottle 1 liter bottle	PowerLyzer™ Tube Holder Stand	13157	1 unit
KAPA PROBE FAST qPCR Kits	51220-100 51220-500 51220-1000	100 reactions 500 reactions 1000 reactions	PowerVac™ Mini System	11992	1 unit + 20 adapters
KAPA SYBR® FAST Universal 2X qPCR Master Mix	51230-100 51230-500 51230-1000	100 reactions 500 reactions 1000 reactions	PowerVac™ Manifold	11991	1 unit
KAPA2G Robust HotStart ReadyMix	51240-100 51240-500	100 reactions 500 reactions	PowerVac™ Mini Spin Filter Adapters	11992-10 11992-20	10 adapters 20 adapters
KAPA HiFi HotStart ReadyMix	51250-100 51250-500	100 reactions 500 reactions	Ceramic Bead Tubes, 1.4 mm	13113-50	50 bead tubes
KAPA2G FAST HotStart DNA Polymerase with dNTPs	51260-100 51260-250 51260-500	100 reactions 250 reactions 500 reactions	Ceramic Bead Tubes, 2.8 mm	13114-50	50 bead tubes
KAPA2G FAST HotStart ReadyMix	51270-100 51270-500	100 reactions 500 reactions	Glass Bead Tubes, 0.5 mm	13116-50	50 bead tubes
KAPA Long Range HotStart DNA Polymerase with dNTPs	51280-100 51280-250 51280-500	100 reactions 250 reactions 500 reactions	Glass Bead Tubes, 0.1 mm	13118-50	50 bead tubes
KAPA Taq Polymerase ReadyMix	51290-250	250 reactions	Metal Bead Tubes, 2.38 mm	13117-50	50 bead tubes
OmniTaq™ DNA Polymerase Enzyme	1224-250	250 reactions (10 U/μl)	2.0 ml Tough Tubes with Cap	13119-500 13119-1000	500 1000
OmniTaq™ DNA Polymerase 2x Master Mix	1226-250	250 reactions (5 x 1.25 ml/tube)	Carbide Bead Tubes, 0.25 mm	13121-50	50 x 0.5 ml tubes
Omni KlenTaq™ DNA Polymerase Enzyme	1225-250	250 reactions (25 U/μl)	Garnet Bead Tubes, 0.15 mm	13122-50	50 x 0.5 ml tubes
Omni KlenTaq™ DNA Polymerase 2x Master Mix	1227-250	250 reactions (5 x 1.25 ml/tube)	Garnet Bead Tubes, 0.70 mm	13123-50	50 x 2 ml tubes
DNase (RNase-Free)	15600-5 15601-100	5 mg 2500 units	Garnet + ¼ Ceramic 15 ml Bead Tubes, 0.70 mm	13134-50	50 tubes
Proteinase K	1223-100 1222-2	100 mg 2 ml (20 mg/ml)	Garnet + ¼ Ceramic 50 ml Bead Tubes, 0.70 mm	13144-10 13144-50 13144-100 13144-500	10 tubes 50 tubes 100 tubes 500 tubes
Ribonuclease A (25 mg/ml)	1202-1 1202-5	1 ml 5 ml	Glass 15 ml Bead Tubes, 0.1 mm	13135-50	50 tubes
PCR Water	17000-1 17000-5 17000-10 17000-11	1 ml 5 x 1 ml 10 x 1 ml 10 ml bottle	Glass 50 ml Bead Tubes, 0.1 mm	13145-10 13145-50 13145-100 13145-500	10 tubes 50 tubes 100 tubes 500 tubes
Molecular Biology Grade Water	17012-200 17012-5200	200 ml 5 x 200 ml	Glass 15 ml Bead Tubes, 1.0 mm	13136-50	50 tubes
DEPC Treated Water	17011-200 17011-5200	200 ml 5 x 200 ml	Ceramic 15 ml Bead Tubes, 1.4 mm	13137-50	50 tubes
Endotoxin-Free Water	17013-10 17013-50 17013-100 17013-500	10 ml 50 ml 100 ml 500 ml	Ceramic 50 ml Bead Tubes, 1.4 mm	13147-10 13147-50	10 tubes 50 tubes
			Metal 50 ml Bead Tubes, 2.38 mm	13149-10 13149-50	10 tubes 50 tubes



## Other Quality Products Available from MO BIO Laboratories, Inc.

For more product and detailed information go to [www.mobio.com/catalog-request](http://www.mobio.com/catalog-request) to request a catalog.

Instrumentation and Accessories... Continued	Catalog No.	Quantity
PowerMix 15 ml Bead Tubes	13138-50	50 tubes
PowerMix 50 ml Bead Tubes	13148-10 13148-50	10 tubes 50 tubes
2 ml Collection Tubes	1200-100-T 1200-150-T 1200-250-T	100 tubes 150 tubes 250 tubes
2 ml Screw Cap Tubes	12800-200-E	200 tubes & caps
15 ml Collection Tubes	12700-T	25 tubes
50 ml Centrifuge Tubes	12600-T	25 tubes
Spin Filters (in 1.9 ml tubes)	1200-50-SF 1200-100-SF 1200-250-SF	50 filters 100 filters 250 filters
Endotoxin-Free Centrifuge Tubes	12617-100  12618-50  12619-25	100 each/2 ml tubes  50 each/15 ml tubes  25 each/50 ml tubes
15 ml Midi Spin Filters	12700-SF	25 spin filters
Vortex-Genie® 2 Vortex (120V)	13111-V	1 unit
Vortex-Genie® 2 Vortex (220V)	13111-V-220	1 unit
Vortex Adapter, holds 12 (1.5-2.0 ml) tubes	13000-V1	1 unit
Vortex Adapter, holds 6 (5 ml) tubes	13000-V1 -5	1 unit
Vortex Adapter, holds 4 (15 ml) tubes	13000-V1 -15	1 unit
Vortex Adapter, holds 2 (50 ml) tubes	13000-V1 -50	1 unit
Vortex Adapter, holds 24 (1.5-2.0 ml) tubes	13000-V1 -24	1 unit
BagMixer® 400 VW	23112	1 unit
BagFilter® 400 P	23113-500	Box of 500
BagPage® 400	23114-500	Box of 500

Instrumentation and Accessories... Continued	Catalog No.	Quantity
Whirl-Pak® Collection Bag, Medium (1,627 ml)	23211-500	500 bags
Whirl-Pak® Collection Bag, Large (3,637 ml)	23212-250	250 bags
Whirl-Pak® Stand up Bag, Small (118 ml)	23220-500	500 bags
Whirl-Pak® Stand up Bag, Medium (532 ml)	23221-500	500 bags
Whirl-Pak® Stand up Bag, Large (1,242 ml)	23222-250	250 bags
Whirl-Pak® Stand up Bag, Extra-Large (2,041 ml)	23223-250	250 bags
Whirl-Pak® Scoop Bag, 60 ml	23240-50	50 bags
Anti-Static Funnels, Micro	23301-96	Pack of 96
Anti-Static Funnels, Small	23302-50	Pack of 50
Anti-Static Funnels, Medium	23303-50	Pack of 50
Anti-Static Funnels, Large	23304-20	Pack of 20
Mini Horizontal Gel System	16001	1 each
Mini Horizontal Gel Caster, 3 place	16003	1 each
Mini Horizontal Gel Tray	16004	1 each
Polycarbonate Single-sided Comb	16005 16006 16007 16008	1 mm x 3 well 1 mm x 8 well 1 mm x 10 well 1 mm x 12 well
Polycarbonate Dual-sided Comb	16013  16014  16015  16016	1 mm x 8 well/16 well 1 mm x 10 well/14 well 2 mm x 8 well/16 well 2 mm x 10 well/14 well
Teflon Single-sided Comb	16009 16010 16011 16012	1 mm x 3 well 1 mm x 8 well 1 mm x 10 well 1 mm x 12 well
Teflon Dual-sided Comb	16017  16018  16019  16020	1 mm x 8 well/16 well 1 mm x 10 well/14 well 2 mm x 8 well/16 well 2 mm x 10 well/14 well
Power Supply w/Timer, (120V)	16023	1 unit



## Other Quality Products Available from MO BIO Laboratories, Inc.

For more product and detailed information go to [www.mobio.com/catalog-request](http://www.mobio.com/catalog-request) to request a catalog.

Instrumentation and Accessories... <i>Continued</i>	Catalog No.	Quantity
96 Well Plate Shaker (120V)	11996	1 unit
96 Well Plate Shaker (220V)	11996-220	1 unit
Plate Adapter Set	11999	1 set

Instrumentation and Accessories... <i>Continued</i>	Catalog No.	Quantity
Vacuum Pump (120V)	11998	1 unit
Vacuum Pump (220V)	11998-220	1 unit